



KOALA
Medical Pty Ltd

EST. 1993

H100B Pulse Oximeter

Adult finger probe or Neonate wrap




Compact Pulse Oximeter ideal for spot checks and transportation. Comes Complete with Adult Probe and Silicone Protective Boot.

The H100B is also Pole Mountable when using the Mount Kit.

Neonatal Mode with Separate Wrap Style Probe Available.

- LCD Display
- Pulse Tone
- Audible and Visible Alarm
- Long Battery Life up-to 50 hours using 4 AA batteries
- Backlight Control Automatic Power Off for Power Saving
- SpO₂ and Pulse Rate Measurement
- Display Numeric Value and Pleth Waveform Simultaneously
- Data Storage Capacity which can be downloaded to a PC
- 10 Minute Review of SpO₂ and Pulse Rate Display Function
- Compatible with Nellcor (OXIMAX) and BCI DB9 sensors

ITEM	CODE
Pulse Oximeter with Silicone Boot and Adult Probe	ED H100B
Adult Pulse Oximeter Probe	ED-H100B Adult Probe
Neonate Wrap Pulse Oximeter Probe	ED-H100B Neonate Wrap
Pole Mounting Kit	ED-H100B Mount Kit



02 9979 9374



sales@koalamedical.com.au



koalamedical.com.au

