




“Maxi-Lite” Shielded Single Use Blades

“Maxi-Lite” Shielded single use laryngoscope blades have been packed with features to ensure there is no compromise in performance.



- Medical grade stainless steel with green spec ISO fitment.
- Shielded light guide, giving maximum light where it is needed most.
- Strong plastic heel, designed to stop the blade tip from touching the handle, eliminating a source of cross-infection.
- Unique twin packaging that can be used as a handle cover, reducing cross-infection.

MACINTOSH CURVED	CODE	QUANTITY
No. 2 Child Macintosh	SG-LL7112LW	Box of 10
No. 3 Adult Macintosh	SG-LL7113LW	Box of 10
No. 4 Large Adult Macintosh	SG-LL7114LW	Box of 10
No. 5 X-Large Adult Macintosh	SG-LL7115LW	Box of 10



02 9979 9374



sales@koalamedical.com.au



koalamedical.com.au

