




**KOALA**  
Medical Pty Ltd

EST. 1993

## Septum Co-Axial Circuit (Centre Divided Single Limb Circuit)



- A versatile single limb circuit capable that reduces clutter and plastic waste, suitable for both adult and paediatric patients.
- These circuits weigh up to 25% less than your regular circle circuit, reducing the torque on the ETT tube and airways.
- Heat exchanged between limbs provides warmer inspired gas.
- Available in both 1.8m and 3.0m lengths.

PRODUCT	DESCRIPTION	CODE	QUANTITY
Septum Circuit 1.8m	Coaxial septum circuit with 22mm divided hose approximately 1.8m with blue extendable rebreathing limb and 2.0 litre rebreathing bag. Patient connector, inspiratory T-piece + adapters.	UD-KO-46658	Box of 20
Septum Circuit 3.0m	Coaxial septum circuit with 22mm divided hose approximately 3.0m with blue extendable rebreathing limb and 2.0 litre rebreathing bag. Patient connector, inspiratory T-piece + adapters.	UD-KO-46657	Box of 20



02 9979 9374



sales@koalamedical.com.au



koalamedical.com.au

