



Sigmoidoscope



The Evexar Self-Illuminating Sigmoidoscope is a single use, fully disposable device that utilises a convenient (universal) Illuminator inserted into the base of the handle, and a light rod to direct light to the end of the scope where it is most needed.

This enables examination and therapeutic procedures to be performed in clinics and the operating room, giving complete clinical freedom of use.

Sigmoidoscopes and Illuminators are completely disposable as a unit, after use.

- Ideal for Specialist Clinicians for use in Hospitals and Private Clinics, as well as for General Practitioners in their surgeries
- Self-illuminated; no need for fibre optic cables or cleaning
- Integrated obturator design prevents obstruction of the examination channel and light guide
- Fully disposable as a unit after use
- Reduced risk of cross-infection since the scope is single use and disposable
- Convenient and cost effective

CODE	DESCRIPTION	ORDER QUANTITY
EV-01111	Sigmoidoscope with Illuminator	Box of 25



DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409



Email: sales@koalamedical.com.au



Proctoscope 19mm and 25mm

The Evexar Self-Illuminating 19mm and 25mm Proctoscopes are single use, fully disposable devices that utilise a convenient (universal) illuminator inserted into the base of the handle, and a light rod to direct light to the end of the scope where it is most needed.

This enables examinations and minor therapeutic procedures to be performed in all areas such as Out-Patient Departments, Accident and Emergency, General Practitioner surgeries, as well as in the operating room, giving complete clinical freedom to use anywhere.



Proctoscopes and Illuminators are completely disposable as a unit, after use.

- Ideal for General Practitioner Surgeries and Out - Patient Clinics; wherever examinations are performed
- Ideal for Specialist Clinicians for use in Hospitals and Private Clinics
- Self-illuminated; no need for fibre optic cables or cleaning
- Integrated obturator design prevents obstruction of the examination channel and light guide
- Fully disposable as a unit after use
- Reduced risk of cross-infection since the scope is single use and disposable
- Convenient and cost effective

CODE	DESCRIPTION	ORDER QUANTITY
EV-01107	19mm Proctoscope with Illuminator	Box of 25
EV-01112	25mm Proctoscope with Illuminator	Box of 25

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au

