



Single Use Connectors



GM-KO-73017 and GM-KO-73018 (with port)

- Single use swivel wye feature, 2 x 22mm
- Limb connection and 22/15mm patient
- Connector ISO fittings and are available with or without Luer Port



UD-KO-72211 and UD-KO-72015 (with port)

- Single use mask elbow 15mm
- Male and female 15mm connections
- Connector ISO fittings and are available with or without Luer Port



GM-KO-76228

- Single use T-Piece with luer port 15mm
- Male and 22mm/15mm male ISO fittings



UD-KO-1423

- Single use metered dose elbow
- Male and 22mm/15mm male ISO fittings

PRODUCT	CODE	QUANTITY
Swivel Wye with Luer Port	GM-KO-73018	Box of 50
Swivel Wye	GM-KO-73017	Box of 50
Mask Elbow with Luer Port	UD-KO-72015	Box of 100
Mask Elbow	UD-KO-72211	Box of 100
Tee-Piece with Luer Port	GM-KO-76228	Box of 50
Metered Dose Elbow	UD-KO-1423	Box of 100

