



KOALA
Medical Pty Ltd

EST. 1993

Therm-O-Finder



NON-CONTACT IR THERMOMETER

- Reading from temporal artery.
- Quick and accurate non-contact temperature detection.
- No covers required with no risk of virus transmission.
- Back-lit display for night use.
- Auto power OFF.

SPECIFICATIONS

- Accuracy +/- 0.2 C
- Operating temperature 15-40 C
- Measurement distance 2-3cm from forehead centre.
- Body temperature measurement range 34.0-42.5 C
- Response time = 1 second.



HOLD 2CM TO
3CM FROM
TEMPLE AREA

PRODUCT	CODE	QUANTITY
Therm-O-Finder non-contact thermometer	FS-300	1 Each



02 9979 9374



sales@koalamedical.com.au



koalamedical.com.au

