



Video Laryngoscope with Reusable Blades



HM-RMAC Blades



HM-VL3R

- Ergonomic lightweight design, anti-microbial handle.
- High resolution video function with internal 4G SD Memory Card, able to take photos and record videos for training/education
- 3.5 inch full colour, non-glare display, convenient for users to observe the laryngeal structure
- 2M pixel excellent Optics Camera, anti-fog design with six LED light, CMOS camera chip
- Micro USB port – recharging or file uploading
- 316 medical stainless steel reusable blade, easy for sterilization by hydrogen peroxide gas plasma or cleaner, 1000 times
- Various size blades available

ITEM	CODE	ORDER QUANTITY
Reusable Video Laryngoscope with 3 Reusable Blades (MAC 2, MAC 3, MAC 4)	HM-VL3R	Each
Reusable Blade MAC 1	HM-RMAC1	Each
Reusable Blade MAC 2	HM-RMAC2	Each
Reusable Blade MAC 3	HM-RMAC3	Each
Reusable Blade MAC 4	HM-RMAC4	Each
Reusable Blade MAC D5	HM-RMACD5	Each
Reusable Blade MILLER 00	HM-RMIL00	Each
Reusable Blade MILLER 0	HM-RMIL0	Each
Reusable Blade MILLER 1	HM-RMIL1	Each

DISTRIBUTED BY

Koala Medical Pty Ltd

Phone: 02 9979 9374

Fax: 02 9979 6409

Email: sales@koalamedical.com.au

