



KOALA
Medical Pty Ltd

EST. 1993

Video Laryngoscope - Disposable Blades



- Ergonomic lightweight design, anti-microbial handle.
- High resolution video function with internal 4G SD Memory Card, able to take photos and record videos for training/education.
- 3.5 inch full colour, non-glare display, convenient for users to observe the laryngeal structure.
- 2M pixel excellent optics camera, anti-fog design with six LED lights with CMOS camera chip.
- Micro USB port – for recharging or file uploading.

- 316 medical stainless steel reusable blade, easy for sterilization by hydrogen peroxide gas plasma or cleaner 1,000 times.
- Various size blades available.

PRODUCT	CODE	QUANTITY
Video Laryngoscope with 3 Disposable Blade packs (MAC2, MAC3, MAC4)	HM-VL3RD DII	1 Set
Disposable Blade MAC 1	HM-DMAC1 DII	1 Each
Disposable Blade MAC 2	HM-DMAC2 DII	1 Each
Disposable Blade MAC 3	HM-DMAC3 DII	1 Each
Disposable Blade MAC 4	HM-DMAC4 DII	1 Each
Disposable Blade MAC D5	HM-DMAC5 DII	1 Each

 02 9979 9374

 sales@koalamedical.com.au  koalamedical.com.au

