



Vista Plus

Transilluminator - Easy Identification Of Veins



The Vista Plus is a portable vein illumination device for both paediatric and adult patients, and patients with a history of difficult venous access.

Vista Plus also improves the visualization of veins in geriatric patients and patients with range of pigmented skin types.

When activated it reveals a map of the vasculature on the surface of the skin. This leads to a more effective venepuncture procedure, with less patient discomfort, resulting in higher patient satisfaction and ensuring shorter procedure time.

- Easy to use, economical and hands free
- Increases the success rate of first attempt venepuncture
- 3 illumination settings for maximum patient versatility
- A rechargeable lithium battery provides 40 min of continuous run time at full illumination & recharges in 2hrs.
- Special indicators show the remaining charge.
- Splash proof, so that the user can clean it with disinfectant

Did You Know?

- More than 50% of children admitted to a healthcare facility requires the insertion of a peripheral intravenous cannula (PIVC).
- More than 50% of these insertions fail on the first attempt.

SIZES AVAILABLE

When ordering please quote the following (minimum order qty 1 box):

ORDER#	DIMENSIONS	QTY/box
TRANS2	162mm x 40mm x 35mm	1

Conlett

U25/17 Cairns St
Loganholme
QLD 4129
Australia

Tel: (+61) 7 3209 9957
Email: contact@conlett.com