




KOALA
Medical Pty Ltd

EST. 1993

Powered Atomization Sprayer



EZ-Spray, the powered atomization sprayer

Koala's EZ-spray was designed for rapid anesthetization of patients. It is a single-patient, closed system that helps to prevent the possibility of cross-contamination. Adjust the oxygen source to operate 7 to 10 LPM according to desired flow rate. The 180' directional tip provides the ability to direct and deliver medication to the upper pharynx effortlessly. It is a quick and thorough anesthetization process.

PRODUCT	CODE	QUANTITY
60cc EZ-sprayer	EZ-100C	Box of 10



02 9979 9374



sales@koalamedical.com.au



koalamedical.com.au

