

PLIANT +

2 Way Single Use Laryngeal Airway



- Complete with compliant and bendable stylet allowing easy insertion & placement
- Made of Medical Grade Silicone, non-toxic & non-stimulating
- Steel reinforced airway tube prevents kinking and occlusion from biting
- Double Layer Cuff seals the throat perfectly
- Gastric Drainage Tube

Reinforced with Gastric Tube- Size 1	<5kg	TM-KO-P1004P2	Box of 10
Reinforced with Gastric Tube- Size 2	10-20Kg	TM-KO-P2010P2	Box of 10
Reinforced with Gastric Tube- Size 3	30-50Kg	TM-KO-P3020P2	Box of 10
Reinforced with Gastric Tube- Size 4	50-70Kg	TM-KO-P4030P2	Box of 10
Reinforced with Gastric Tube- Size 5	>70Kg	TM-KO-P5040P2	Box of 5

“Right Here, Right Now”



KOALA
Medical Pty Ltd